




# High-throughput pH Measurement

Olin College SCOPE for moderna


# Lipid nanoparticle formulations deliver mRNA to the body

- Each mRNA product requires a different LNP formulation
- There's no model to predict performance
- Experimental validation is necessary




# Moderna needs to characterize thousands of LNP formulations

*pH is one of the most basic electrochemical measurements and is relevant to both the stability and biocompatibility of drug products*



# Single probes are rate-limiting


Because no multichannel probes exist, even robotic solutions are still very slow (~90 minutes to measure 96 samples).






# Screen Printed Electrodes (SPEs) offer benefits over glass probes


- Fluids can be pipetted directly onto the surface
- Many sensors in one device
- No risk of crashing




pH measurement requires a selective element to distinguish between ions



# No ion-selective electrode card existed



# Now one does!








We built a pH sensor that can take 96 measurements in 15 minutes.



# System Overview




# Retaining System




**Retaining System**  
**Electrode Card**  
**Interface Board**  
**Measurement Board**  
**Software**

# Electrode Card



Retaining System  
**Electrode Card**  
Interface Board  
Measurement Board  
Software


# Interface Board



Retaining System  
Electrode Card  
**Interface Board**  
Measurement Board  
Software



# Measurement Board



Retaining System  
Electrode Card  
Interface Board  
**Measurement Board**  
Software

# Software

```
# show -p
```

```
Showing pH values from the most recent measurement in pH units. To see the associated voltages, use 'show -v'.
```

	1	2	3	4	5	6	7	8	9	10	11	12
A	10.00	10.00	10.00	4.00	10.00	10.00	10.00	10.00	4.00	10.00	10.00	4.00
B	4.00	10.00	10.00	4.00	10.00	10.00	10.00	4.00	10.00	10.00	4.00	10.00
C	9.81	10.00	8.12	4.00	4.00	10.00	5.18	10.00	10.00	4.00	10.00	4.00
D	4.00	10.00	10.00	10.00	5.23	4.00	4.00	10.00	10.00	5.38	10.00	9.16
E	10.00	5.50	4.00	10.00	10.00	10.00	4.00	10.00	10.00	4.00	10.00	10.00
F	4.00	4.62	10.00	10.00	10.00	4.00	4.00	10.00	4.00	4.00	4.00	10.00
G	4.89	10.00	4.00	6.81	10.00	10.00	10.00	5.79	10.00	4.00	10.00	10.00
H	4.00	4.00	4.00	4.00	10.00	4.00	4.15	10.00	4.00	7.39	4.00	10.00

```
# |
```


**Retaining System  
Electrode Card  
Interface Board  
Measurement Board  
Software**

# Our Device

- Increases throughput by 30x
- Works with microvolume samples
- Removes manual data entry
- Enables use of a liquid handling robot

## and allows Moderna to...

- Characterize more LNP solutions
- Ultimately develop mRNA medicines faster



# Thank You!

- Thank you to our wonderful advisor, Alessandra
- Thank you to folks at Moderna: Rud, Xavier, Stephanie, Nick, Andrew

